


How to confirm the most suitable irradiation time?

1. At your first treatment, choose the irradiation time to 30 seconds, if your skin doesn't appear skin erythema or feel the burn, increase 20 seconds of the irradiation time on the second day.
2. Follow the above process, and each irradiation time increased 20 seconds until the skin appears pink mild erythema.
3. Record the treatment time, then reduce 20 seconds of that, that's your most suitable irradiation time.
4. After you make sure your best treatment time, using once every 2-3 days, the treatment effect will be better.

Notice:

During the treatment, if time above 5 minutes, your skin still doesn't appear erythema or feel the burn, we suggest do not increase the irradiation time, the longest irradiation time is 5 minutes once.

Operation instruction



① "Mode"option : manual and automatic :

Manual mode : After starting, the timer is counting along with irradiation,

user needs to turn off the lamp by yourself;

Automatic mode : After starting, the timer is countdown as your setting. (auto turn-off)

② "Time"option :

MM stands minute, SS stands second.

③ "Start"option :

Press the button for 3-4 seconds to turn on or turn off the lamp,
when you hear "beep" sound, release the button.

Special function :

After starting, no option is selected in the case:

Press the modify button ⑦ 3 seconds to switch the language.

Press the modify button ⑤ 3 seconds to switch left hand to right hand.

④ ⑤ ⑥ ⑦ "Select" button :

④ ⑥ "Select" button, press the button to select a different option, the

item is selected when blue background light stopped.

⑤ ⑦ "Modify"button, press this button to modify the contents of the blue background light option. Such as adjusting the exposure time, ⑦ Increase

the number, ⑤ Reduce the number.