

Product Manual of Beauty Equipment

Quick Start Guide:

- Put on the laser protective glasses before starting any laser treatment.
- Ensure the glasses fit securely and comfortably on your face.
- Adjust the elasticity to accommodate different face shapes.
- Remove the glasses after the treatment is complete.

Warnings:

- CHOKING HAZARD - Small parts. Not for children under 3 yrs.
- Consult a doctor before use if you have any eye conditions.

Safety Information:

- Keep dry and store in a cool, dry place.
- Wipe with a damp cloth to clean the lenses.

Service: Contact our store support for any issues.