

USER MANUAL

1. Product Info

Product Name:Laser Machine

Model: TM001

Batch Number

2. MANUFACTURER & EU REPRESENTATIVE INFO

- Manufacturer: Changsha Haoyu International Trade Co., Ltd, Room 1641, Building 1, 9, 2, Zone 2, Jinqiao Market Cluster, South of Pu Ruixi Road, Huangjinyuan Street, Xiangjiang New
- District, Hunan, yanglulu.lu99@gmail.com/+86 18890551572
- EU Authorized Representative: EU Consulten Srls, Via Varesina 55 Milano 20156 Italia Campobasso
- Molise,Eu-consultenSrls@hotmail.com/39-3285832654

3. Safety & Warnings

- *Protective eye goggles must be worn during operation.*
- *Contraindicated groups: pregnant women, people with photosensitive disorders, open skin wounds and epilepsy sufferers.*
- *This device is only designed for hair removal on arms, legs, underarms and bikini areas. Do not use on scalp, nipples. Facial use is prohibited unless certified for facial treatment.*
- *Ensure cooling vents are unobstructed during use to avoid overheating.*
- *Do not operate near flammable gas.*
- *CE mark, EU authorized representative information, batch number and model are permanently printed for compliance check.* Warning: Children, Keep the device and its cord out of reach of children under 3 years old.

4. Instructions for Use

如：

1.Preparation: Shave hair thoroughly on target skin areas; wipe skin dry with no skincare products left. Remove all protective films on the laser window and body.

2.Power Supply: Plug the matched power adapter into standard household socket and connect to the device.

本模版由速卖通提供，仅供参考。卖家应根据产品的实际特性及销往国的法律要求自行调整。速卖通不对因使用本模版产生的任何合规或法律问题承担责任。This template is provided by AliExpress for reference only. Sellers should adjust the content based on actual product features and local laws. AliExpress assumes no responsibility for any compliance or legal issues arising from the use of this template.

3. *Power On:* Press and hold the power button for 3 seconds to boot up the machine. Put on laser safety goggles before any operation.

4. *Modes & Energy:* Short press the function button to switch different energy levels. Start with low energy for first use.

5. *Operation:* Fit the laser window tightly against flat skin; the device emits laser automatically upon complete contact. Slide the device slowly for full coverage.

6. *Power Off:* Long press the power button for 3 seconds to shut down. Unplug the adapter after the device cools down.

5. MAINTENANCE & STORAGE

如:

- *Cleaning:* Turn off the device before cleaning. Use a soft, lint-free cloth.
- *Storage:* Store in a dry place at room temperature (15°C - 35°C).
- *Long-term Storage:* If the device will not be used for a long time, fully power off and unplug. Clean the lens thoroughly. Avoid squeezing the power cord sharply..

6. DISPOSAL & RECYCLING



This symbol indicates that this product should not be treated as household waste. According to the WEEE Directive (2012/19/EU), it must be disposed of at an appropriate electronic waste collection point for recycling to prevent potential negative consequences for the environment and human health.

7. TECHNICAL SPECIFICATIONS

Feature	Laser / Intense Pulsed Light
---------	------------------------------

本模版由速卖通提供，仅供参考。卖家应根据产品的实际特性及销往国的法律要求自行调整。速卖通不对因使用本模版产生的任何合规或法律问题承担责任。This template is provided by AliExpress for reference only. Sellers should adjust the content based on actual product features and local laws. AliExpress assumes no responsibility for any compliance or legal issues arising from the use of this template.

Feature	Laser / Intense Pulsed Light
Input Voltage	110-220V ~ 50Hz
Rated Power	10W
Battery Capacity	2000mAh
Certifications	CE, RoHS, UKCA

本模版由速卖通提供，仅供参考。卖家应根据产品的实际特性及销往国的法律要求自行调整。速卖通不对因使用本模版产生的任何合规或法律问题承担责任。This template is provided by AliExpress for reference only. Sellers should adjust the content based on actual product features and local laws. AliExpress assumes no responsibility for any compliance or legal issues arising from the use of this template.