


195x100mm
2025.11.26




MULTI FUNCTIONAL
CLAMPING MASSAGER

QUICK START GUIDE

MULTI FUNCTIONAL
CLAMPING MASSAGER

Disclaimer: Any self-care plan carries inherent risks.
Before starting this or any treatment plan,
please consult with your healthcare professional for physical therapy.
his product is not intended to treat any medical condition.
We are not responsible for any harm or damage caused by use.
Use at your own risk. This product is not intended for any other use.

GET TO KNOW
YOUR MASSAGER



1.ROLL CAGE (MASSAGE ROLLERS)
2.POWER INDICATOR LIGHTS
3.POWER BUTTON

4.RED LIGHT SWITCH
5.CHARGING PORT
6.HAND GRIPS

HOW TO USE

Be sure to read all warnings & instructions before operating



Red Light Switch
Pushing the red light switch upwards turns the red light on and pushing it downwards turns it off.


Power On
Firmly press and hold the power button on the top of the handle for about 2 second. Roll cage will start spinning automatically.

Set Speed
Each short press of the power button sets the speed (1-8) as indicated by the lights above.

Reverse scrolling
Double click the power button twice to scroll in the opposite direction

Power Off
Firmly press and hold the power button on the top of the handle until the device turns off.

COMPLETE BODYWELLNESS



BEFORE USING OR CHARGING THE MASSAGER, READ ALL THE INSTRUCTIONS, WARNINGS, AND CAUTIONARY MARKINGS IN THIS MANUAL.

The massager is intended for home use only and not by children.
Do not use the massager on your head, face, or near your genitals.
Consult your doctor before using this if you are pregnant.
Do not allow loose hair or clothing near any moving part of the massager.
Do not use the massager above your throat or vertebrae.
Do not get the massager wet. Only clean the massager with a lightly damp cloth.
Do not cover the device with anything and do not block the vents of the motor.
Never attempt to disassemble the massager.
Never store the massager at temperatures above 40°C/104°F. Only charge the massager when the ambient temperature is between 0°C/32°F and 40°C/104°F. If the device feels hot during use, immediately cease use.
Never use or charge the massager near water or expose to water.
For long-term storage, the massager should be fully charged, but not stored plugged into the charger.
This product is not intended to diagnose, treat, or prevent any diseases, In the unlikely event you begin to experience discomfort or pain during use, immediately cease use and if needed, consult a doctor.
While the massager is unlikely to cause injury, there is a small likelihood that it could exacerbate existing conditions or injuries. You assume all risk of use of the massager.

Warning: Please Tie Up Your Long Hair Before Use!hair Entraming May Cause Serious Injury!