Instructions

002

004

Operation and use method

1.002 How to turn on

- ※ Plug in the Android mobile phone data cable. After power on, the instrument enters the standby state, and the light indicator is always on.
- ※ Long press the "switch" key on the instrument for 1 second to start the machine, the instrument starts to work, the indicator is always on, and turn on conversion vibration massage.
- ※ The intensity of conversion vibration massage is divided into three grades, which will be automatically adjusted to the highest grade after startup. After startup, briefly press the "switch" key to adjust the gear.
- ※ Press and hold the "switch" key on the instrument for 1 second to turn it off.

2.002 How to connect the power supply

- ※ Please use the USB charging cable provided with the instrument to connect the mobile phone charger or charging treasure (the output voltage is 5V and the current is not greater than 1.5A). After power on, the light is on and the light is on after work.

3.002 Why doesn't the product have a charging adapter

- ※ The product supports universal Android USB data cable
- ※ Support the standard universal mobile phone adapter and powerbank on the market.
(Note: Output Voltage 5V, output current > 5A, recommended output current > = 1a)

1.004 How to turn on

- ※ Plug in the Android mobile phone data cable. After power on, the instrument enters the standby state, and the light indicator is always on.
- ※ Press the "switch" key on the instrument to start up, the instrument starts to work, the indicator light is always on, starts vibration massage, left and right massage, conversion vibration. ※ Press the "switch" key on the instrument to shut down.
- ※ During use, press the "+" key once to increase the strength by one gear, Each time you press the "-" key, the strength will be reduced by one gear. There are 3 gears in total. Please choose your own comfort intensity.

2.004 How to connect the power supply

- ※ Please use the USB charging cable provided with the instrument to connect the mobile phone charger or charging treasure (the output voltage is 5V and the current is not greater than 1.0A). After power on, the light of the indicator is on and the light is on after work.

3.004 Why doesn't the product have a charging adapter

- ※ The product supports universal Android USB data cable
- ※ Support the standard universal mobile phone adapter and powerbank on the market.
(Note: Output Voltage 5V, output current > 5A, recommended output current > = 1a)

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