

W Air Compression Leg Massager

User Guide

WE'RE HERETO HELP!

ThankriGils
e
described in this User Guide may void the limited warranty.
Please do not touch the unit.
This is it.

Box Contents
• Leg massager
• Adapter
• User guide

Parts
1. Leg massager
2. Air pressure hose
3. LCD Screen
4. Control panel



Safety Precautions

- Do not use voltage beyond the requirement of this product.
- Do not attempt to repair or dismantle the product. There are no user-serviceable parts.
- Do not use the product in a place where it may become wet or damp. If the unit becomes wet, do not use.
- Stop using the product immediately if you feel uncomfortable or any abnormalities occur.
- Do not use this product in a high-temperature environment.
- Do not use this product if you are currently pregnant or menstruating.
- Do not place the unit in a place where it may be in contact with flammable items, explosive items, or chemicals.
- Do not touch the unit with any sharp object or risk damage.
- Do not use this product if any of the below apply to you:
 - Skin is infected, swollen, burnt, or otherwise damaged.
 - You have a pacemaker or other medical implant.
 - You are suffering from heart disease.
 - Your bones are fragile or injured.
 - You are suffering from cancer.
 - You are receiving treatment or your blood pressure is not normal.
 - You are suffering from osteoporosis.

WARNING:
Use only as instructed. Do not disassemble or modify the product. Keep away from water and fire sources.

FUNCTION DESCRIPTION OF FACE BUTTONS



HAND CONTROL OPERATION INSTRUCTION BOOK

- Plug in the power adapter (the power hole is the lowest CD socket hole), all the function buttons are white;
- Click the main unit power-on switch 1 and all the function buttons are blue, the main unit enters the intelligent massage mode, and then click the main unit to turn off the main unit and turn white.
- Click the massage mode selection button 2, the corresponding 1, 2, and 3 modes are displayed in the indicator box 13, and the 3 modes click cycle;
- 4, 3, and 4 are the keys for adjusting the pressure of massage strength, the strength is 3 levels, which can only be adjusted in one direction, and the corresponding level is displayed in the 13 number box during adjustment;
- Click the massage part selection button 5 to switch between feet, calves, and full legs, and display correspondingly in the indicator areas 10, 11, and 12, and 11 and 12 are at the same time for full legs. Display, click the massage mode selection button 2 to exit the massage part selection and enter the full-leg mode;
- Click the hot compress function selection button 6, click once to enter the low-grade L, click again to enter the high-grade H, and click again to turn off the hot compress function, the corresponding file The position is correspondingly displayed in the gear display area 13.

Operation

- Slip on each massager, keeping your feet flat on the bottom. Note that each is marked L and R for each foot.
- Use the strap to tightly secure your foot into the massager.
- Use the strap to tightly secure your leg into the massager.
- Move your foot to make sure that the massager is securely fastened.
- Insert the two air pressure connectors into the controller.
- Press the power button. While on the LED will stay on.
- Adjust the air pressure strength and mode by your hand, and can see your operation in LCD screen.

NOTE: the device will automatically shut off after 20 minutes.
For best results, do not use the leg massager for more than an hour.

Care and maintenance

- Store the product in a dry place, away from moisture and direct sunlight.
- Wipe clean with a dry towel.

Troubleshooting

Problem	Cause	Solution
The device does not work.	The unit is not connected to a power source.	Make sure the unit is fully connected to the power.
The power button was not pressed.	The power button was not pressed.	Press the power button.
The device suddenly stopped working.	The unit was running more than 15 minutes.	Restart the unit to resume use.

Specifications

Rated voltage: 12V 2A
Power consumption: 24W
Default time: 15 min