

# PRODUCT MANUAL

## 1. Product Informaton

**Product Name:** [Teeth whitening lamp]

**Model:** [764]

**Batch/Serial Number:** [05]

## 2. MANUFACTURER & EU REPRESENTATIVE INFO

- **Manufacturer:** [SHENZHEN CHOICURATION FUTURE TECHNOLOGY CO., LTD.]
- **Address:** CN, 205, Floor 2, No.1, Chenghe Industrial Zone, Tongsheng Community, Dalang Street, Longhua District, Shenzhen City, Guangdong Province, 44031;
- **Email:** 1604299951@qq.com;
- **Tel:** 86-13537565865]
- **EU Authorized Representative:** [Authorized Representative Company];
- **Name:** APEX CE SPECIALISTS LIMITED;
- **Full Address:** UNIT 3D NORTH POINT HOUSE, NORTH POINT BUSINESS PARK, NEW MALLOW ROAD, CORK, T23AT2P, LRELAND;
- **Email:** info@apex-ce.com;
- **Tel:** 353-212066339]

## 3. Compliance Information



本模版由速卖通提供，仅供参考。卖家应根据产品的实际特性及销往国的法律要求自行调整，以确保最终用户手册**准确、完整并符合法律要求**。速卖通不对因使用本模版产生的任何合规或法律问题承担责任。This template is provided by AliExpress for reference only. Sellers should adjust the content based on actual product features and local laws to ensure that the final user manual is **accurate, complete, and legally compliant**. AliExpress assumes no responsibility for any compliance or legal issues arising from the use of this template.

*NOTE: This product must comply with all applicable safety and regulatory requirements in the market where it is sold. Sellers must assess which requirements and regulations apply based on the actual product type, function, and destination market.*

#### **4. Safety & Warnings**

- This product is a home teeth whitening aid and is only suitable for healthy adults. It is contraindicated for minors, pregnant women, and breastfeeding women. It is also contraindicated for those with sensitive teeth, periodontal disease, enamel damage, or oral lesions. Consult a dentist before use.
- Do not look directly at the LED light source to avoid blue light damage to the eyes. Wear protective goggles or close your eyes while using this product.
- This product is for personal use only and must not be shared with others to avoid cross-infection. Do not disassemble, modify the lamp body, or replace non-original parts.
- Do not immerse the device in water or charge/use it in a humid environment. Use only the included data cable when charging to avoid damaging the circuitry.

#### **5. Instructions for Use [Mandatory]**

1. Preparation before use: Clean and dry the teeth whitening tray; connect the teeth whitening lamp to the mobile phone/USB power source and confirm that the lamp is emitting light normally.
2. Instructions for Use: Apply the whitening gel evenly to the inside of the mouthpiece as instructed in the product instructions. Wear the mouthpiece and turn on the device. Recommended usage time is 10-15 minutes per session; avoid exceeding this time.
3. Cleaning and Maintenance: After use, turn off the power, remove the mouthpiece, rinse thoroughly with water, and dry. Store in a cool, dry place. Clean the lamp surface regularly to prevent oil and dust from obstructing the light source.

*NOTE: Sellers shall add and adjust product-specific steps where necessary.*

本模版由速卖通提供，仅供参考。卖家应根据产品的实际特性及销往国的法律要求自行调整，以确保最终用户手册**准确、完整并符合法律要求**。速卖通不对因使用本模版产生的任何合规或法律问题承担责任。This template is provided by AliExpress for reference only. Sellers should adjust the content based on actual product features and local laws to ensure that the final user manual is **accurate, complete, and legally compliant**. AliExpress assumes no responsibility for any compliance or legal issues arising from the use of this template.