

2 in 1 Diode Laser Hair Removal Machine 1064 755 940

Thank you for choosing our product. Please read the following instructions carefully before use.

Product Name: Multi-Function Personal Care Device

Model: Standard

Origin: Mainland China

Power Source: USB Cable / Rechargeable Battery / Dry Battery / Electric / Dry & Rechargeable / USB Charger / Outside Charger / Wireless Charger

Voltage: 100V / 110V / 220V / 240V (Please confirm your local voltage before use)

Batteries Included: Yes

Material: High-quality ABS and electronic components

Working Principle: Normal Light / ION+ / ION- / Low Level Laser / LED / RF / EMS / Warm / Cold / Infrared / Ultrasound Wave / Vibration Massage

Functions: Lifting / Whitening / Blemish Removal / Skin Moisture / Anti-Acne / Skin Rejuvenation / Facial Clean / Anti Wrinkle / Skin Tightening / Dead Skin Removal

Certifications: CE / RoHS / FCC / UL / GS
/ ETL / PSE / KC / PSB / CCC

Place of Origin: Mainland China

Power Plug Options: US / EU / UK / AU /
Japanese / Korean

Medical Supervision Required: No (For
personal and home use only)

Hand Instrument: No

Usage Instructions:

1. Charge the device fully before first use. Use the included USB cable or charger.
2. Clean the treatment area and ensure skin is dry and free of oils or lotions.
3. Select the desired function and intensity level using the control button.
4. Gently move the device over the target area in a circular motion.
5. Do not use on broken or irritated skin.
6. Recommended usage time: 10-15 minutes per session, up to 3 times per week.

Safety Precautions:

- Keep away from water and high humidity environments.
- Do not disassemble or modify the device.
- Stop use immediately if skin irritation or discomfort occurs.
- Keep out of reach of children.
- Do not expose to extreme temperatures or direct sunlight.

Maintenance:

- Clean the device head with a soft, dry cloth after each use.
- Store in a cool, dry place.
- Recharge the battery at least once every three months if not in use.

For any questions or support, please contact our customer service team.